




neodisher® TP alka

designed for Tecniplast



Detergent for the automated reprocessing of animal cages

Liquid concentrate

Fields of application:

- Automated cleaning of animal cages and accessories as well as cage, storage and transport racks in specially designed washers

Performance Spectrum:

- Alkaline, intensive detergent with a high proportion of dispersants
- Reliably removes organic residues such as faeces, blood and proteins
- Compatible with stainless steel, polycarbonate, polyetherimide, polysulfone and polyphenylsulfone

Special Properties:

- Free of phosphates and surfactants
- Without oxidizing agents
- Excellent material compatibility

Application and dosage:

neodisher TP alka is used in specially designed washers for animal cages and accessories as well as cage, storage and transport racks. The dosing amount also depends on the degree of contamination of the items to be washed and is 2 - 6 ml/l.

Use suitable dosing devices for dosing.

When using neodisher TP alka the water for the cleaning step should not exceed a total water hardness of 0.5 mmol/l.

General instructions on use:

- The neodisher TP alka working solution has to be rinsed off completely (preferably with deionised water).
- Do not mix with other products
- Rinse out dosing system including suction hose with water before changing product.
- The instructions given by the manufacturer of the washers are to be observed.
- Please observe the recommendations on reprocessing given by the manufacturer of the items to be washed and the recommendations of the working group for cage processing” (AK KAB) in the current issue of the AK KAB brochure “Cage Processing in Animal Facilities properly done”.
- For professional use only

Technical data:

pH-range	12.4 - 12.9 (2 - 6 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C):
Density	1.4 g/cm ³ (20 °C)
Titration factor	0.37 (in accordance with neodisher titration method)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:
5 - 15 % polycarboxylates

Storage information:

Always store at a temperature ranging between 0 and 30 °C.

Usable for 4 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hour glass symbol



neodisher[®] TP alka

designed for Tecniplast

Hazard warnings and safety advice:

Classification and labelling according to the
Dangerous Preparation Directive
1999/45/EC for the concentrate as delivered.

Hazard symbol C - corrosive

contains: potassium hydroxide

Risk phrases:

R 22 - Harmful if swallowed.

R 35 - Causes severe burns.

Safety phrases:

S 26 - In case of contact with eyes, rinse
immediately with plenty of water and seek
medical advice.

S 27/28 - After contact with skin, take off
immediately all contaminated clothing, and wash
immediately with plenty of water.

S 36/37/39 - Wear suitable protective clothing,
gloves and eye/face protection.

S 45 - In case of accident or if you feel unwell,
seek medical advice immediately (show the label
where possible).

Dispose only when container is empty and closed. For
disposal of product residues, refer to Material Safety
Data Sheet.

For further safety information see EC-Safety data
sheets. These are available e.g. at www.drweigert.de
under the headline "Service".

Sales units	Contents
Drum	280 kg
Canister	25 kg
Canister	10 L

MB 4082/3-1 Revision Date 06/2012

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.