




neodisher® TP oxy

designed for Tecniplast



Component for the automated cleaning of aquatic tanks

Liquid concentrate

Field of application:

- Automated cleaning of aquatic tanks in specially designed washers

Performance spectrum:

- Additional component based on active oxygen
- Reliably removes stubborn residues such as algae, food residues and water ingredients in combination with neodisher TP alka
- Suitable for tanks made of polycarbonate and polysulfone

Special properties:

- Free of active chlorine and surfactants
- Easy to rinse off
- Excellent material compatibility

Application and dosage:

neodisher TP oxy is used in combination with the alkaline detergent neodisher TP alka in special cleaning units for aquatic tanks. Among other things the dosage depends on the degree of soiling of the aquatic tanks and is 3 - 5 ml/l.

Use suitable dosing devices.

When using neodisher TP oxy the water for the cleaning step should not exceed the total hardness of 3 °d (0.5 mmol/l).

General instructions on use:

- The neodisher TP oxy working solution has to be rinsed off completely (preferably with deionised water).

- Do not mix with other products.
- Rinse out dosing system including hose with water before changing product.
- For professional use only.
- The instructions given by the manufacturers of the cleaning units are to be observed.
- Please observe the recommendations on reprocessing given by the manufacturers of the aquatic systems.
- For professional use only.

Technical data:

pH-value	5.4 – 5.3 (3 - 5 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	1.1 g/cm ³ (20 °C)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:
 > 30 % oxygen-based bleaching agents

Storage information:

Always store at a temperature between -20 °C and + 25 °C. Do not expose to direct sunlight.

Usable for 1 year when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol ⏳.



neodisher[®] TP oxy

designed for Tecniplast

Hazard and precautionary statements:

For safety and environmental information see EC safety data sheets.

Dispose only when container is empty and closed.
For disposal of product residues, refer to Material Safety Data Sheet.

MB Stand: 4085/3-2

Revisions date: 04/2015

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.