






# neodisher® TP alka

designed for Tecniplast



## Detergent for the automated reprocessing of animal cages

Liquid concentrate

### Fields of application:

- Automated cleaning of animal cages and accessories as well as cage, storage and transport racks in specially designed washers
- Automated cleaning of aquatic tanks in specially designed washers

### Performance Spectrum:

- Alkaline, intensive detergent with a high proportion of dispersants
- Reliably removes organic residues such as faeces, blood and proteins
- Compatible with stainless steel, polycarbonate, polyetherimide, polysulfone and polyphenylsulfone

### Special Properties:

- Free of phosphates and surfactants
- Without oxidising agents
- Excellent material compatibility

### Application and dosage:

neodisher TP alka is used in specially designed washers for animal cages and accessories as well as cage, storage and transport racks. The dosing amount also depends on the degree of contamination of the items to be washed and is 1 - 6 ml/l.

Cleaning fish tanks requires specially matched processes which makes contacting the manufacturer of the specially designed washers necessary in order to determine the correct dosage.

Use suitable dosing devices.

When using neodisher TP alka the water for the cleaning step should not exceed a total water hardness of 0.5 mmol/l.

### General instructions on use:

- The neodisher TP alka working solution has to be rinsed off completely (preferably with deionised water)
- Do not mix with other products
- Rinse out dosing system including suction hose with water before changing product
- The instructions given by the manufacturer of the washers are to be observed
- Please observe the recommendations on reprocessing given by the manufacturer of the items to be washed and the recommendations of the working group for cage processing" (AK KAB) in the current issue of the AK KAB brochure "Cage Processing in Animal Facilities properly done"
- For professional use only

### Technical data:

pH-range	12.0 - 12.8 (1 - 6 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C):
Density	1.4 g/cm <sup>3</sup> (20 °C)
Titration factor	0.35 (in accordance with neodisher titration method)

### Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:  
5 - 15 % polycarboxylates




# neodisher<sup>®</sup> TP alka

## designed for Tecniplast

### Storage information:

Always store at a temperature between 0 °C and 30 °C.

Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

### Hazard and precautionary statements:

For safety information see EC safety data sheets.

Dispose only when container is empty and closed.  
For disposal of product residues, refer to the Safety Data Sheet.

MB 4082/3-3    Revision Date 11/2018



---

With the above information, to our current knowledge we describe our product regarding safety necessities, but we do not involve any quality description or promise certain properties.

---