




neodisher[®] TP rinse

designed for Tecniplast



Rinsing agent for the automated reprocessing of animal cages



Fields of application:

- Rinsing of animal cages and accessories as well as cage, storage and transport racks in specially designed washers

Performance Spectrum:

- Considerably improves and accelerates drying
- Suitable for polysulfone, polycarbonate, polyphenylsulfone und polyetherimide

Special Properties:

- Excellent material compatibility, prevents the formation of stress cracks during autoclaving
- Toxicologically safe
- Acidic, neutralises alkaline residues

Application and dosage:

neodisher TP rinse is used in specially designed washers for animal cages and accessories as well as cage, storage and transport racks. The dosing amount also depends on the material of the items to be washed and is 0.5 – 2.0 ml/l.

Use suitable dosing devices.

For the automated cleaning of animal cages and accessories we recommend using neodisher TP alka in the alkaline process and neodisher TP acid in the acidic process.

General instructions on use:

- Do not mix with other products.
- Rinse out dosing system including hose with water before changing product.
- For professional use only.
- The instructions given by the manufacturer of the washer are to be observed.
- Please observe the recommendations on reprocessing given by the manufacturers of the animal cages and the items to be washed as well as the recommendations of the working group for cage processing (AK KAB) in the current issue of the AK KAB brochure "Cage Processing in Animal Facilities properly done".

Technical data:

pH-range	approx. 3.6 – 3.2 (0.5 - 2 ml/l, determined in deionised water, 20 °C)
Density	approx. 1.04 g/cm ³ (20 °C)
Viscosity	< 10 mPa s (concentrate, 20 °C)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:
< 5 % polycarboxylates
also preservatives
(chloromethylisothiazolinone/
methylisothiazolinone)

Storage information:

Always store at a temperature between 0 °C and 30 °C.
Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol ⌚.

Hazard and precautionary statements:

Dispose only when container is empty and closed.
For disposal of product residues, refer to Safety Data Sheet.

For further safety information see EC safety data sheets.

MB 4084/3-3 Revision Date 11/2019