

**TECNIPLAST**  
*innovation through passion*



**Sealsafe**

IVC BLUE LINE



HOUSING



THE NEW **"LIFT & LOCK" SYSTEM** IS DESIGNED TO HELP OPERATORS WHEN CHANGING A CAGE. THANKS TO THE NEW LID DESIGN, THE CAGE TOP CAN BE LOCKED IN THE OPEN POSITION TO REACH ANIMALS MORE EASILY. MOREOVER, CAGES REQUIRE SIGNIFICANT LESS SPACE ALLOWING THE 2000P IVCS TO FIT UNDER MOST STANDARD CHANGING STATIONS AND LAMINAR FLOW EQUIPMENT.

- Great visibility-focus design: easy daily health check may be carried out with no need to move cages and disturb animals
- Ergonomics: cages can be pulled out partially, in case of closer inspections
- Latches, gasket and microbiological filter guarantee a high level of protection for operators and animals

### VERSATILITY

- The divider allows to split 2000P cage in two halves to run different experiments and enhance density
- The fully accessible internal cage height, together with the cage optional mezzanine, permit animals to play, jump and rest, providing a more comfortable three dimensional environment compared to other standard caging systems
- Versatile system for housing either rats or guinea pigs

### ANIMAL WELFARE

- Air inlet and outlet valves are located in the top: no air draft at animal level, avoiding risk of stress and heat loss

### HIGH PROTECTION FOR OPERATORS AND ANIMALS

- Nylon frame with silicone co-moulded gasket on the lid enhances the sealing between cage and filter top
- Microbiological filters (VFE $\geq$ 99.999987% and BFE  $\geq$ 99.9999937%) at cage level, ensure optimal gas exchange as well as protection during procedures or in case of power failure



LIFT & LOCK SYSTEM



EASY AND QUICK BOTTLE CHANGE/FOOD REFILLING



PARTIAL CAGE SLIDING OUT

Cage	Rack Configuration
2000P	15 cages (3w x 5h), Single sided

COMPANY WITH  
QUALITY SYSTEM  
CERTIFIED BY DNV GL  
= ISO 9001 =

COMPANY WITH  
ENVIRONMENTAL SYSTEM  
CERTIFIED BY DNV GL  
= ISO 14001 =

**TECNIPLAST**  
innovation through passion

#### HEADQUARTERS

ITALY | Tel. +39 0332 809711 · www.tecniplast.it · E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it  
SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. + 44 (0)845 0504556 · www.tecniplast.it · E-mail: info@tecniplastuk.com

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CAT171A/ING - Rev Oct2018

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.