



INNOVATION CAN BE SMALL IN SIZE BUT GREAT IN PRINCIPLE

AUTOCLAVABLE FULL PLASTIC CONSTRUCTION WITH CRUSHPROOF STRUCTURE

The BC-Plast bottle crates are **molded injected in a single piece** using a glass-filled polyphthalamide (PPA) which exceeds standards in heat, moisture and chemical resistance. The unique plastic design, by nature, is guarantee of no bending and proper shape even under intense usage and handling.

PATENTED QUICK LOCK & RELEASE EMBEDDED LID

A patented solution allows to avoid the annoying removable lids used to hold in place the bottles when manipulated upside down. An **embedded grid** permits to **lock or release bottles** with a simple sliding movement; this improves any activity from cage changing to washing especially when combined with the drawer style transport trolleys.

35% LIGHTER THAN EVER


With less than 3 kg (6.5 lbs) of weight, the BC-Plast is at least 35% lighter than any previous S/S crate. That, combined with the easy mechanism to access the bottles, is synonymous with **top in class for ergonomics**.

FLEXIBLE TO ACCOMMODATE DIFFERENT BOTTLE SIZES

The BC-Plast comes as a standard in two possible configurations capable to accommodate bottles from 200ml (13oz) to 400ml (26oz). Please consult us for special custom needs.

SUITABLE FOR MANUAL AND AUTOMATED PROCESSING

This next generation accessory has been developed to be compatible with the most common equipment, from the more manual emptying washing and filling lines to the fully automatized systems. **The BC-Plast can be the constant in your flows of today and tomorrow.**



CLOSE



OPEN



a **TECNIPLAST** company



HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · www.iwtsrl.it · E-mail: info@iwtsrl.it

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. + 44 (0)845 0504556 · www.tecniplast.it · E-mail: info@tecniplastuk.com

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

TO FIND YOUR LOCAL DISTRIBUTOR, PLEASE VISIT WWW.IWTSRL.IT

CAT W-DP065ING Rev May 2016

IWT reserves the right to change or modify product and/or specifications without notice or obligation.