



- Versatile system for housing either guinea pigs or rats
- Improved animal welfare thanks to the limited intra-cage air speed
- Great visibility-focus design: easy daily check
- Latches, gasket and microbiological filter guarantee a high level of protection for operators and animals

## VERSATILITY

- The two available drop-in hoppers allow the use of 2000P IVC either for Guinea Pigs or Rats

## ANIMAL WELFARE

- Air inlet and outlet valves are located in the top: no air draft at animal level, avoiding risk of stress and heat loss
- The fully accessible internal cage volume makes for a more three-dimensional and thus stimulating housing solution, compared to other standard caging systems

## HIGH PROTECTION FOR OPERATORS AND ANIMALS

- Nylon frame with silicone co-moulded gasket on the lid enhances the sealing between cage and filter top
- Microbiological filters (VFE\* $\geq$ 99.999987% and BFE\* $\geq$ 99.9999937%) at cage level, ensure optimal gas exchange as well as protection during procedures or in case of power failure

Cage	Rack Configuration	
2000P	2P15	15 cages (3w x 5h), Single sided



COMPLETE CAGE



LID WITH OPEN HATCH



LID

### HEADQUARTERS

ITALY | Tel. +39 0332 809711 · [www.tecniplast.it](http://www.tecniplast.it) · E-mail: [tecnicom@tecniplast.it](mailto:tecnicom@tecniplast.it); [venditeitalia@tecniplast.it](mailto:venditeitalia@tecniplast.it)

### SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · [www.tecniplast.it](http://www.tecniplast.it) · E-mail: [info@tecniplast.com.au](mailto:info@tecniplast.com.au)

CHINA | Tel. + 86 (0)21 50810920 · [www.tecniplast.cn](http://www.tecniplast.cn) · E-mail: [tecniplastchina@tecniplast.it](mailto:tecniplastchina@tecniplast.it)

FRANCE | Tél. +33 (0)4 72 52 94 41 · [www.tecniplast.fr](http://www.tecniplast.fr) · E-mail: [info@tecniplast.fr](mailto:info@tecniplast.fr)

GERMANY | Tel. +49 (0)8805 921320 · [www.tecniplast.de](http://www.tecniplast.de) · E-mail: [info@tecniplast.de](mailto:info@tecniplast.de)

UNITED KINGDOM | Tel. + 44 (0)845 0504556 · [www.tecniplast.it](http://www.tecniplast.it) · E-mail: [info@tecniplastuk.com](mailto:info@tecniplastuk.com)

JAPAN | Tel. +81 (0)3 5770 5375 · [www.tecniplast.it](http://www.tecniplast.it) · E-mail: [info@tecniplastjapan.co.jp](mailto:info@tecniplastjapan.co.jp)

USA | Toll Free: +1 877.669.2243 · [www.tecniplastusa.com](http://www.tecniplastusa.com) · E-mail: [info@tecniplastusa.com](mailto:info@tecniplastusa.com)

CANADA | Toll Free: +1 855.347.8718 · [www.tecniplastcanada.com](http://www.tecniplastcanada.com) · E-mail: [info@tecniplastcanada.com](mailto:info@tecniplastcanada.com)

CAT170/ING - Rev Dec2016

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.