



Effective algae and biofilm removal
No need for scrubbing or tank pretreatment

HIGH EFFICIENCY and COMPACT DESIGN

FITS YOUR NEEDS. FITS IN YOUR LAB.

- No need for manual scrubbing or tank pretreatment
- Greater detergent savings and reduced running costs
- Simplified utilities - No need for exhaust duct
- "Zero gravity" vertical door to washing chamber
- Modular building for easy rigging and installation

PRODUCT FEATURES

- Suitable for tanks and accessories
- Vertical door at ergonomic height
- Wide see-through glass door for full load inspection
- POLARIS interface for remote connectivity
- Multiple customizable washing cycles
- Externally validated disinfection process
- Final rinse with reverse osmosis (RO) water
- Frontal access to technical chamber
- On-board compartment to store detergents

- Maximum Exterior Installed Dimensions (frontal view)
53" W x 35" D x 78" H
(1350mm x 892mm x 1980mm)
Minimum recommended ceiling height 2300mm

- Maximum Rigging Dimensions (frontal view)
53" W x 33" D x 77" H
(1350mm x 835mm x 1950mm)



BENEFITS

- **Utilizes 80% less** detergent than prior models
- **Simplified** install - No need for exhaust duct
- **Compact** dimensions - Fits standard door width
- **Reduced** operational footprint
- **Ergonomic** height for tank loading and unloading
- **Wall-to-wall installation** - Frontal access to maintenance chamber



Independently validated software-driven process delivers perfectly clean and residue-free tanks