

*The most reliable and versatile  
single use caging solution made  
of recycled plastic designed  
to satisfy specific needs in terms of:*

- **Transportation**
- **Quarantine**
- **Contingency and Disaster planning**
- **Temporary Solutions**
- **Irradiation Studies**
- **Satellite Facilities**
- **Company Start up/Spin off**



## 1/The Solution for mice transportation

The confirmed scientific importance and success of transgenic mice has led to thousands of different colonies being shipped from one research lab to another. SUMC represents the solution for animal transportation: once at the animal receiving area it can be docked to the IVC Tecniplast Sealsafe Plus standard rack by mean of simple adaptors for housing or quarantine purposes.

## 2/Flexibility: one cage for multiple purposes

Tecniplast further increases its wide range of product solutions with the SUMC system, a 500 square centimeter disposable housing solution made of recycled plastic that provides the traditional ventilated cage performance and TP reliability you've come to expect from us.

SUMC can be the ideal solution for specific needs such as:

- Quarantine
- Stocking contingency caging in case of an emergency
- Temporary housing solutions
- Irradiation studies
- Satellite facilities with limited washing capabilities
- Work in biotech start-up companies with limited resources

The same cage for static applications may be used in transportation and IVC housing, thanks to the adaptors.

The high stackability makes it easy to store and manage.

Our system offers the following benefits:

Features	Benefits
Transparent Plastic Material: A-PET	Environmentally friendly, maximizing end use recycling. Better visibility
Stackable	Easy storage, nesting densely with 115 cage bases in 1 meter, reducing storage space
One product for multiple options: transportation/static/IVC	High Flexibility: designed to work in static as well as an IVC in existing Tecniplast racks, adopting the superior TP method of air delivery at top level
Versatile new water bottle	Better condition for mice also during transportation

SUMC components are available sterile (Gamma Irradiated) double bagged or not sterile single bagged.

## Cage Liner

Tecniplast single use product range is completed with Polyethylene cage liners.




# HOUSING

## SUMC - Single Use Mouse Cage

Flexibility

1&2



# TECNIPLAST

HEADQUARTERS

ITALY | Tel. +39 0332 809711 • www.tecniplast.it • E-mail: tecnicom@tecniplast.it

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 • www.tecniplast.it • E-mail: info@tecniplast.com.au

FRANCE | Tel. 04 72 52 94 41 • www.tecniplast.fr • E-mail: info@tecniplast.fr

GERMANY | Tel. 08805 921320 • www.tecniplast.de • E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. 0845 0504556 • www.tecniplast.it • E-mail: info@tecniplastuk.com

JAPAN | Tel. + 03 5770 5375 • www.tecniplast.it • E-mail: info@tecniplastjapan.co.jp

USA/CANADA | Toll Free: 877.669.2243 • www.tecniplastusa.com • E-mail: info@tecniplastusa.com

To find your local distributor, please visit [www.tecniplast.it](http://www.tecniplast.it)



ISO-9001 - Cert. n° 0875



ISO-14001 - Cert. n° 0256A/0



Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation. SUMC is a Tecniplast trademark.