











## **Top-of-the-range Ventilated Cabinet**

Tecniplast ARIA BIO-C36 EVO is a Ventilated Cabinet 1.0 m (36") in size, designed to accommodate nearly all types of cages with no cross contamination.

## PERFORMANCE

- Innovative dual ventilation system with two energy saving blowers in order to work in positive or negative pressure according to the users' needs
- Microprocessor driven ACH with an automatic compensation of the air flow rate, to keep safe working conditions over time
- Supply and Exhaust air filter HEPA H14 efficiency (EN 1822) with DOP test scanning design.

## BENEFITS

- Flexibility: Positive (for bioexclusion) and Negative (for biocontainment pressure mode settings
- Cleanability: exhaust air ducts and shelves simple to dismount and clean
- Energy-efficient and Acoustic Comfort due to the unique plenums design and EC motors
- Optional innovative and dimmable LED lighting system, with a daily/ weekly timer, throughout the whole surface of each compartment
- Optional reliable heating system with uniform and efficient heat transfer.

## STANDARD FEATURES

- User friendly 7" TFT Colour Touch panel with integrated microprocessor system for easy control and operation
- Dual-paned magnetic tempered safety glass doors available in transparent, red or opaque color
- Working area entirely built in AISI 304 S/S
- Adjustable shelving for flexible cage housing (only for cabinets without day/night cycle option)
- Compliant with CE requirements and UL/CSA.





OVERALL DIMENSIONS (WxDxH) mm 1000 x 810 x 1930 29" 3/8 x 31" 7/8 x 76"

WORKING AREA - without shelves- (WxDxH) mm 940 x 680 x 1350 37" x 26" 6/8 x 53" 1/8



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001=

COMPANY WITH
ENVIRONMENTAL SYSTEM
CERTIFIED BY DNV GL
= ISO 14001=



**HEADQUARTERS** 

ITALY | Tel. +39 0332 809711 · www.tecniplast.it · E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com SWITZERLAND | Tel. +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com CAT197/ING · Rev June2021

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.