



UP & DOWN





SPACE-SAVING AND FLEXIBLE DESIGN

- **Light design** and **compact dimensions** for minimal footprint and high manoeuvrability
- **Plug&Play** system with no need of compressed air
- **Electrical lifting system** for bedding reservoir:
 - low noise solution
 - no bedding material quality loss

DISPENSING EFFICIENCY AND VERSATILITY

- Up to **350 Typell cages** for each silo refill
- Two independent dispensing points with horizontal augers for **high dispensing accuracy**
- Up to **10 different pre-settable volumes** for each point
- Full flexibility to process **different cage sizes**
- A **wide range of bedding types** and a set of selected **environmental enrichments** can be processed

OPERATOR SAFETY AND ERGONOMICS

- **Working height fully adjustable**
- **Suction ring** technology for high-level operator protection during refilling and dispensing operations
- **Double level filters (G4-F7)** and efficient airflow for dust and small particles containment
- **Automatic anti-collision system** during the lowering in resting position



ADJUSTABLE WORKING HEIGHT



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · www.iwtsrl.it · E-mail: info@iwtsrl.it

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

SWITZERLAND | Tel +41 (0)81 750 31 40 · www.indulab.ch · Email: infos@indulab.ch

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CAT MKTBRC000003H Rev March2020

IWT reserves the right to change or modify product and/or specifications without notice or obligation.