



The new Tecniplast portable DeconTank is a **compact solution** specifically designed for **ISOcage decontamination**. The DeconTank is dedicated to all customers that are working with a **Biosafety Cabinet** without an integrated dunk tank, but wish to achieve a **complete and safe decontamination** of the cage external surface.

The DeconTank is a **plastic box with a holding chamber** that allows complete coverage of disinfectant on the cage surface, allowing for **appropriate external decontamination**. All this just by operating a front lever.



## EASY AND SAFE

- **Effectiveness of the decontamination:** no more «home made» solutions that are not guaranteeing a suitable and complete decontamination.
- **Ergonomics: Handy, portable solution** that can be placed on any surface/table to be located near the Biosafety Cabinet. Thanks to the optional trolley, the height is fully adjustable to be aligned to the opening of the BSC. Tap for drain to ease the decontaminant emptying drained.
- Designed to **minimize the amount of disinfectant** used and **avoid spillages**.
- **Cleanability:** the holding chamber can be dismantled for cleaning the inside of the chamber.
- **Control of the contact time** thanks to the integrated timer.

COMPANY WITH  
QUALITY SYSTEM  
CERTIFIED BY DNV GL  
= ISO 9001 =

COMPANY WITH  
QUALITY SYSTEM  
CERTIFIED BY DNV GL  
= ISO 9001 =

**TECNIPLAST**  
innovation through passion

### HEADQUARTERS

ITALY | Tel. +39 0332 809711 · www.tecniplast.it · E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it

### SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tel. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplast.it · E-mail: info@tecniplastuk.com

SWITZERLAND | Tel +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CAT198/ING - Rev Apr2021

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.