

 **TECNIPLAST**



ARIA ISO-TECH60 N
BIOSAFETY CABINET



ARIA

ISO-TECH60 N

Class II Biosafety Cabinet for High-Containment Laboratories

Advanced Safety & Containment Features

- **Dual motor blower system** – A **second motor blower** in the exhaust box ensures containment even in case of a primary blower failure.
- **Airtight front window (optional inflatable gasket)** – Enhances safety and simplifies H₂O₂ in loop decontamination.

Uncompromising Performance & Compliance

- **ISO Class 3 Air Quality** (ISO 14644-1)
- Aeraulic performance in accordance with **EN 12469:2000** as **Class II biosafety cabinet** and **NSF49** as a **Class II Type A2 biosafety cabinet**

Key Features & Benefits

- **Automatic sliding front window** – Two modes for flexibility:
 - Biosafety mode (250 mm front opening)
 - Changing mode (330 mm front opening)
- **7" TFT color touch panel** – Real-time airflow monitoring with an intuitive interface
- **Energy-efficient LED lighting** – Includes dimmable white and red LED options for night operations
- **ECM motor blower with automatic airflow control** – Ensures consistent laminar airflow and air barrier performance
- **TRESPA® worktop panels** – Highly durable and easy to clean, making routine cleaning and autoclaving more efficient

Second motor blower

Manual Decon Lock

Optional inflatable gasket

Automatic Dunk Tank

TECNIPLAST

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

COMPANY WITH
ENVIRONMENTAL SYSTEM
CERTIFIED BY DNV GL
= ISO 14001 =

HEADQUARTERS

ITALY | Tel. +39 0332 809711 · www.tecniplast.it · E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.jp · E-mail: info@tecniplastjapan.co.jp

SWITZERLAND | Tel +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com