






## The Tecniplast Isolator is designed for facilities working with Germ Free or Axenic animals, that wish to breed or maintain stock animals while preserving the experimental part in ISOCages.

Isolator provides bioexclusion for gnotobiotics, germ-free or immunocompromised animals due to its design to maintain a sterile environment with a constant supply of high positive pressure, clean HEPA









## SAFETY & PROTECTION

- Supply and exhaust air is pre-filtered and HEPA filtered guaranteeing clean air
- High positive pressure inside the Isolator canopy
- Transfer chamber equipped with a double air lock rounded port to transfer supplies or connect to other equipment to maximize the protection
- Pressure alarm to notify if the pressure has decreased under a certain threshold
- Dry-contact to connect alarms to the Building Monitoring System.

## DURABILITY

- High-quality flexible canopy wrapped on a stainless steel frame to avoid any risk of collapse and damage resulting from power failures
- Removable sleeves with easy change telescopic cuffs: both the gloves and sleeves can be individually changed
- Stainless steel main structure provides strength and protection for the canopy while moving the isolator.

## ERGONOMIC AND USER FRIENDLY

- Additional service gloves to ease the procedures of assembling and setting up the Isolator
- Two sizes available for different requirements: 2 gloves, for 15 cages or 4 gloves, for 30 cages
- Mini rack specifically designed for ISOCage P cage bases allows easy transportation from Isolators to the ISO rack
- Oval sleeves ensure comfort for different operators' height.

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

COMPANY WITH
ENVIRONMENTAL SYSTEM
CERTIFIED BY DNV GL
= ISO 14001 =



**HEADQUARTERS** 

ITALY | Tel. +39 0332 809711 · www.tecniplast.it · E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplast.it · E-mail: info@tecniplastuk.com SWITZERLAND | Tel. +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CAT193/ING - Rev Dec2020