





## Quiet and compact dumping station optimizing safety and efficiency

Tecniplast DS One bedding disposal station provides **the best-in-class safety**, **ergonomics** and **energy efficiency** for **personnel protection** and **working environment** from aerosol particles and unpleasant odors.

### PERFORMANCE

- Enhanced operator protection against allergens and contaminants thanks to a redesigned and optimized working area;
- Extended HEPA filter life span guaranteed by the pre-filters (ISO Coarse ≥70% as per ISO 16890; G4 as per EN779) clamped with a transparent co-molded net for an easier and simpler removal of larger dust particles;
- Sliding reduction funnel allowing the operator to simplify waste bag removal under a constant and incomparable operator protection.

#### BENEFITS

- Spacious working area in a compact and ergonomic design thanks to its unmatchable access opening height. An optional electric lifting system is also available to be adjusted for user convenience:
- Easiest Cleanability thanks to the working surface made in high tech, highly resistant plastic making it highly cleanable and with a high resistance to chemical disinfectants;
- Energy-efficient and Acoustic Comfort owed to the unique plenum design and EC motor;
- Increased durability of the working surface. The elastic property of the plastic (polymer) is extremely durable and resistant to chemical disinfectants;
- Top-mounted exhaust carbon filter absorbs odours and gaseous contaminants including ammonia

#### STANDARD FEATURES

- AISI 304 with scotch brite finishing for enhanced cleaning and product life;
- 4.3" touch screen control panel with automatic regulation system to keep the inflow air speed constant;
- LED light Bar with light dimmer provide the highest lighting level;
- EC Motor Blower ensures higher efficiency & performance, extended life and reduced maintenance.

# Wide range of options and accessories:

- 60L Polyethylene bin with S/S trolley
- -105L Polyethylene bin with S/S trolley
- Stainless steel waste bag trolley
- Thimble duct exhaust transition
- Exhaust carbon filter
- Trespa shelf





OVERALL DIMENSIONS (WxDxH) with fixed stand support:

mm 1050 x 895 x 1880 - 41" 3/8 x 35" 2/8 x 74" with lifter:

mm 1050 x 895 x 1880 / 2050 - 41" 3/8 x 35" 2/8 x 74" / 80" 6/8

WORKING AREA DIMENSIONS (WxDxH) mm 966 x 660 x 590 - 38" x 26" x 23" 2/8





COMPANY WITH QUALITY SYSTEM CERTIFIED BY DNV GL

COMPANY WITH
ENVIRONMENTAL SYSTEM
CERTIFIED BY DNV GL
= ISO 14001 =



**HEADQUARTERS** 

ITALY | Tel. +39 0332 809711 · www.tecniplast.it · E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp SWITZERLAND | Tel +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CAT187/ING - Rev Apr2020

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.