


The EMber is a system designed to increase temperature inside Emerald IVCs.

It provides a **comfortable, warm environment** to animals for many different application areas such as post surgical recovery, metabolism of drugs, body thermoregulation studies, drug therapy and metastasis and cancer investigations. It can also be used to help weak cage inhabitants, such as pups or specific strains, to maintain homeostasis.

Heating automatically starts when **the cage equipped with the Metal Plate** is placed in the rack, thanks to the **induction heating process** generated by the **Heating Module**. **Temperature** can be set up to 40°C and, once reached, it is **continuously monitored and strictly controlled, keeping it constant**.

EM500 standard IVC racks can be equipped with multiple EMbers, which are available for both standard cages and EM500 featuring the InRichment.




HEATING MODULE



CAGE BODY WITH METAL PLATE

EMber is very easy to use thanks to **the user friendly interface**, which allows the desired experimental conditions to be **intuitively set**.



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

COMPANY WITH
ENVIRONMENTAL SYSTEM
CERTIFIED BY DNV GL
= ISO 14001 =

TECNIPLAST
innovation through passion

HEADQUARTERS
ITALY | Tel. +39 0332 809711 · www.tecniplast.it · E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it
SISTER COMPANIES
AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au
CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it
FRANCE | Tel. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr
GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de
UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplast.it · E-mail: info@tecniplastuk.com
SWITZERLAND | Tel +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch
JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp
USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com
CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CAT206/ING - Rev Mar 2025
Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.