





EOLUS automated Bedding Handling System is designed for the continuous handling, transportation and disposal of dirty bedding materials. EOLUS increases safety and efficiency, and improves workflows.

WASTE DISPOSAL STATION (WDS)

- **Working flexibility:** compatible with all common types of bedding and enrichment available on the market thanks to the shredder technology.
- **High-level of operator protection for minimum allergen exposure:** back draft hood with a combination of G4+F7+H14 filters.
- **Efficient material transportation:** material flow controlled by the conveyor belt which also ensures enrichment processability.



WASTE COLLECTION SYSTEM

- **Maximum flexibility** for different waste collection systems (*four wheeled mobile bins, sealed container and open dumpster solutions*): system customized for your specific workload.
- **Easier and safer operations:** bag clamp, suction ring technology, easy camlock connections (*for sealed container solutions*).
- **Filter lifespan optimization:** embedded pre-filtration system for enrichment and fibrous bedding applications.

TECHNICAL COMPARTMENT

- **Compact dimensions with minimal footprint:** compact, all-in-one solution which includes the vacuum pump(s), the filtration unit and the control box.
- **Layout flexibility:** available in three sizes (COMPACT, REGULAR, DOUBLE) and designed to fit all specific customer requirements.
- **Filtration unit high-level operator protection:** with easy-sealing bag system and HEPA filter (H14).



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · www.iwtsrl.it · E-mail: info@iwtsrl.it

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

SWITZERLAND | Tel +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CAT MKTBRC000041 Rev Sept2020

IWT reserves the right to change or modify product and/or specifications without notice or obligation.