





EOLUS automated Bedding Handling System is designed for the continuous handling, transportation and dispensing of bedding material. EOLUS increases safety and efficiency and improves workflows.

LOADING FUNNEL + SILO(S)

- **Working flexibility:** modular design to easily switch from tote to small bag applications and vice versa
- **High-level operator protection and safety:** suction rings for dust exposure reduction
- **Efficient system workflow:** single or multiple intermediate silo(s) available in different sizes for larger buffer capability.

DISPENSER UNIT

- **Maximum flexibility:** different dispenser models for both manual and automation applications (EASY, FLEXY and FLEXY XL versions)
- **High-level dispensing versatility:** embedded pre- filtering system and mixer for a wide range of bedding types processing, including light-weight bedding materials
- **High-level operator usability and protection:** integrated touchscreen panel and efficient suction rings technology at the dispensing points.

TECHNICAL COMPARTMENT

- **Compact dimensions with minimal footprint:** compact, all-in-one solution which includes the vacuum pump, the filtration unit and the control box
- **Layout flexibility:** available in two sizes (COMPACT, REGULAR) and designed to fit all specific customer requirements.
- **Filtration unit with high-level operator protection:** easy-sealing bag system and HEPA filter (H14).



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · www.iwtsrl.it · E-mail: info@iwtsrl.it

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

SWITZERLAND | Tel +41 (0)81 750 31 40 · www.tecniplast.ch · Email: info@tecniplast.ch

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CAT MKTBRC000107M Rev March2021

IWT reserves the right to change or modify product and/or specifications without notice or obligation.