



**WASHING**

EASY 200 is a pass-through washing system for efficient cleaning of water bottles and bottle caps used in animal facilities. The unique ergonomic design aims to reduce additional effort and fatigue during every repetitive operator task.

Complete with multi-type process tables including bottle emptying, de-capping, re-capping, and refilling. . . . .



## SPACE-SAVING AND FLEXIBLE DESIGN

- Electrically heated bottle washer with **compact dimensions** to keep a full bottle process in a limited space
- Suitable both for **standalone and in-line pass-through applications** featuring double manual sliding interlocked doors
- **Flexible** to be installed in existing layout with an adaptable workflow direction and configuration
- **Full flexibility to be integrated with multiple tables for operators' ergonomic processing**, including semi-automatic solutions for de-capping and re-capping.

## CLEANING EFFICIENCY AND EFFICACY

- **High washing performance** with separated washing and rinsing circuits and a throughput of 36 bottles/cycle (#2 crates/cycle)
- **Dedicated and fixed S/S nozzles system** composed of #36 separated washing and rinsing jets
- **Tempered glass window** installed on the bottle washer side for efficient process phases monitoring.

## ENVIRONMENTALLY- FRIENDLY IMPACT SOLUTION

- **Re-circulating washing cycle** with consequent reduced water and detergent consumptions
- **Heat recovery system to save up to 34% of energy** during a single cleaning cycle
- **Water discharge control system** (temperature and pH control) for leaner integration with customers' main lines.



COMPANY WITH  
QUALITY SYSTEM  
CERTIFIED BY DNV GL  
= ISO 9001 =

HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · [www.iwtsrl.it](http://www.iwtsrl.it) · E-mail: [info@iwtsrl.it](mailto:info@iwtsrl.it)

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · [www.tecniplast.it](http://www.tecniplast.it) · E-mail: [info@tecniplast.com.au](mailto:info@tecniplast.com.au)

CANADA | Toll Free: +1 855.347.8718 · [www.tecniplastcanada.com](http://www.tecniplastcanada.com) · E-mail: [info@tecniplastcanada.com](mailto:info@tecniplastcanada.com)

CHINA | Tel. + 86 (0)21 50810920 · [www.tecniplast.cn](http://www.tecniplast.cn) · E-mail: [tecniplastchina@tecniplast.it](mailto:tecniplastchina@tecniplast.it)

FRANCE | Tél. +33 (0)4 72 52 94 41 · [www.tecniplast.fr](http://www.tecniplast.fr) · E-mail: [info@tecniplast.fr](mailto:info@tecniplast.fr)

GERMANY | Tel. +49 (0)8805 921320 · [www.tecniplast.de](http://www.tecniplast.de) · E-mail: [info@tecniplast.de](mailto:info@tecniplast.de)

JAPAN | Tel. +81 (0)3 5770 5375 · [www.tecniplast.it](http://www.tecniplast.it) · E-mail: [info@tecniplastjapan.co.jp](mailto:info@tecniplastjapan.co.jp)

SWITZERLAND | Tel +41 (0)81 750 31 40 · [www.tecniplast.ch](http://www.tecniplast.ch) · Email: [info@tecniplast.ch](mailto:info@tecniplast.ch)

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · [www.tecniplastuk.com](http://www.tecniplastuk.com) · E-mail: [info@tecniplastuk.com](mailto:info@tecniplastuk.com)

USA | Toll Free: +1 877.669.2243 · [www.tecniplastusa.com](http://www.tecniplastusa.com) · E-mail: [info@tecniplastusa.com](mailto:info@tecniplastusa.com)

CAT MKTBRC000115J Rev June2021

IWT reserves the right to change or modify product and/or specifications without notice or obligation.