



ARCADIA



WASHING



ADVANCE DESIGN: FUNCTIONALITY MEETS OPERATOR NEEDS

- At least **15% footprint reduction** versus comparable previous configurations
- **Ergonomic loading height** of 3' – 900mm
- **Higher loading capacity** given by a tunnel opening of 18" - 460mm
- **Plastics mould-injected conveyer belt** which allows **IVC tops inclined presentation**



BEST PERFORMANCE HAND IN HAND WITH THE LOWEST CONSUMPTIONS

- Up to **40% water consumptions reduction** with respect to previous comparable models
- **Up to 1,500 mice cages/hour**
- Fully validated microbiological results compliant with AK-KAB and AAALAC requirements
- **Heat-recovery system** (optional) for extra energy efficiency

MAXIMIZE USABILITY, MINIMIZE ORDINARY MAINTENANCE TIME

- **Automatic slide-up service doors** for full chambers accessibility
- **Self-cleaning in-line filter** with back-flush feature
- Polaris Interface with remote connectivity
- Full frontal maintenance for **"against the wall" installation** capability



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · www.iwtsrl.it · E-mail: info@iwtsrl.it

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

SWITZERLAND | Tel +41 (0)81 750 31 40 · www.indulab.ch · Email: infos@indulab.ch

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CAT W-DP069ING Rev Oct2018

IWT reserves the right to change or modify product and/or specifications without notice or obligation.