



IWT SKATEBOARD – THE NEW TROLLEY FOR MANUAL OPERATIONS


The IWT skateboard is the new high-capacity transport trolley specifically designed for manual handling of multitype cage bases. A single accessory **to match the full vivarium operational loop**: from animal holding room to the washing area and back through autoclaving.

AUTOCLAVABLE PLASTIC (PSU) CAGE HOLDERS CONSTRUCTION

High-quality and resistant cage holding kits made of mold-injected polysulfone (PSU). Each set of holders is individually engineered to **ensure cage stack stability during transportation** and handling procedures.

FLEXIBLE DESIGN TO ACCOMMODATE DIFFERENT TYPES OF CAGE BASES

Different cage holder models are provided as additional kits to accommodate a wide range of small rodent cages. The flexible design of the trolley allows **up to 3 different cage types to be accommodated**.



SMART STORAGE, "QUICK&EASY" SETTING

Handy Storing slots for having additional sets of cage holders always with you. **Rapid** handwheel knobs guarantee a fast **trolley transformation**.

MODULAR CONFIGURATION

IWT's skateboard could be equipped with additional features to meet all facility requirements: bumpers, s/s push-pull ergonomic handle, **allergen lid for operators' protection** and coupling kit to chain together multiple trolleys for a safe autoclave loading/unloading.

	Emerald	GM500	GR900	GR1800DD	1145T	1285L	1500U
Cage Bases	96	80	30	30	60	60	10



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · www.iwtsrl.it · E-mail: info@iwtsrl.it

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

SWITZERLAND | Tel +41 (0)81 750 31 40 · www.indulab.ch · Email: info@indulab.ch

UNITED KINGDOM | Tel. + 44 (0)345 0504556 · www.tecniplastuk.com · E-mail: info@tecniplastuk.com

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CAT W-DP070ING Rev May2019

IWT reserves the right to change or modify product and/or specifications without notice or obligation.