

OCEANUS



WASHING

OCEANUS



HIGH EFFICIENCY CABINET WASHER

The **Most Powerful Solution** to Operational Efficiency within a **Reduced Footprint**

- Compact design and simplified installation
- More loading capacity and better ergonomics
- Simplified utility requirements
- Less operational footprint... and much more!

BENEFITS

- 25% less footprint than other washers on the market
- 40% more loading capacity than previous versions
- 40% gain of operational footprint
- "Zero-gravity" vertical door for easy operation
- Wide see-through glass door for full load inspection
- Suitable for cages, bottles and accessories
- Optional exhaust conditioning system for duct-free installations



No pull-out trays for better ergonomics. Flexible configuration to fit cages or bottle crates (metal or the new BC-PLAST)



Utilizes the **Atlantis'** proven high-performance washing and rinsing technology.

STANDARD FEATURES

- Separate washing and rinsing circuits
- Full coverage oscillating arm technology
- Polaris interface for remote connectivity
- USB port for data gathering
- Validated microbiological disinfection

- Full front access to electrical panel
- Wall-to-wall installation
- Double loading shelves with flexible configuration
- Fits up to 32 mouse cages

The most powerful washing solution for small and medium lab animal facilities.

The perfect back up for large facilities.