

TOP ADULT FISH WELFARE

Enhanced schooling behaviour of fish with **increased swimming area** and shallow water
Environmental enrichment by means of gravel bottom background

HEALTHY STOCK MAINTENANCE

On-rack easy and efficient spawning thanks to the **dedicated spawning platform**
Breeding events are easily scheduled without disturbing fish


GREATER WATER QUALITY ALSO DURING SPAWNING

Through the **patented siphon technology**
Innovative **water distribution pipe** with micro holes which produce a laminar water flow at the bottom of the tank

FLEXIBLE

Easy to be retrofitted on existing Tecniplast racks

All this in a
single, flexible
and
retrofittable tank!



16 l
volume

HEADQUARTERS

ITALY | Tel. +39 0332 809711 • www.tecniplast.it • E-mail: tecnicom@tecniplast.it; venditeitalia@tecniplast.it

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 • www.tecniplast.it • E-mail: info@tecniplast.com.au

CHINA | Tel. + 86 (0)21 50810920 • www.tecniplast.cn • E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 • www.tecniplast.fr • E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 • www.tecniplast.de • E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. + 44 (0)345 0504556 • www.tecniplastuk.com • E-mail: info@tecniplastuk.com

JAPAN | Tel. +81 (0)3 5770 5375 • www.tecniplast.it • E-mail: info@tecniplastjapan.co.jp

SWITZERLAND | Tel +41 (0)81 750 31 40 • www.tecniplast.ch • Email: info@tecniplast.ch

USA | Toll Free: +1 877.669.2243 • www.tecniplastusa.com • E-mail: info@tecniplastusa.com

CANADA | Toll Free: +1 855.347.8718 • www.tecniplastcanada.com • E-mail: info@tecniplastcanada.com

CAT182/ING - Rev Apr2020

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

COMPANY WITH
ENVIRONMENTAL SYSTEM
CERTIFIED BY DNV GL
= ISO 14001 =

TECNIPLAST
innovation through passion